Tuta GVS2

GVS2 Coverall

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Tuta GVS2

GVS2 Coverall

Protective coverall against biological/chemical agents and cytotoxics

Technical features

GVS2 coverall is protective clothing, CE marked as Personal Protective Equipment (PPE) in category III for the protection against biological and chemical agents, in accordance with requirement of EC Directive 686/89.

This coverall, having ultrasonic welded and taped seams, has obtained the EC type – examination certificate in accordance with the requirements of harmonized standard EN 14126/2004, EN 14605, EN13982-1, EN340, EN13034 for the protection against biological, chemical agents and cytotoxics.

Therefore this PPE is suitable to protect against biological and chemical agents with reference to the EC Directive 679/90 and to EC Directive 54/2000.

GVS2 coverall can also be used for the protection against chemical agents as acids, bases, solvents, disinfectants, cytotoxics, in sanitary and other fields in compliance to the EC Directive 24/98 and EC Directive 394/90. PPE patented by CL.com: WO 03/041807A1

Special Technical features

GVS2 coverall is made up of two different materials, assure a perfect wearing and comfort allowing normal movements, are physiologically innocuous and hydrophobic.

Latex free, manufactured with ultrasonic and taped seams.

Application/areas of use

This GVS2 Coverall can be used by the operator in those working activities when there is exposure or potentially exposure to biological agents as bacteria, parasites, fungi, viruses (biological agents of group 2 and 3 in reference to Annex III of EC Directive 54/2000) as for example HBV, HCV, HIV, and other agents responsible of Creutzfeldt-Jakob, BSE, Avian Flu, etc. and for the protection against chemical agents as acids and bases, as well as for the protection against cytotoxics.

Certification

Protection against biological agents: verified by Notified Body in accordance with EN 14126/04 and in accordance with American Society for Testing Materials ASTM F 1670/97 – F 1671/97

standards.

Protection against chemical agents: verified by Notified Body in accordance with the requirements of harmonized standard EN 14605, EN 13982-1, EN 340 and EN 13034 concerning the tests to be run using chemical agents as strong acids and bases generally used in industrial field and in hospital settings.

The worker protection from exposure to cytotoxics has been checked from the Notified Body, using in the test provided by the above mentioned harmonized standard

Models:

SFT00201- Hooded coverall and with ultrasonic seams reinforced for protection against biological and chemical agents (chemotherapy / cytotoxic)

SFT00301- Hooded coverall and with ultrasonic seams reinforced for protection against biological and chemical agents (chemotherapy / cytotoxic) with shoes and stop-cuff

Tuta GVS2

GVS2 Coverall

Materials (latex free & solvent free)

Polypropylene 40 g/m² + polyethylene 25 g/m²

Color: White

Sterile

The PPE can be sterilized by ethylene oxide

Storage

The correct storage conditions are: room temperature from + 5 °C to + 40 °C ; R.H. < 60 %

The PPE must be stored in a cool and dry place. They cannot stand the sunshine, fluorescent lights, x-rays, ozone and damp.

Expire date

5 years from production date, following the above mentioned storage conditions in the original packaging.

Instruction for use

Please before use read the Instruction for use

Application limits

This clothing does not protect the operator against viral agents in group 4, with reference to the EC Directive 54/2000.

This garment does not protect against gas or areic.

This clothing does not protect against heat/fire, cold, electricity (flammable stuff if getting in touch with fire or sources of heat) and against other kinds of risk not explicitly mentioned.

This clothing only protects the covered parts of the body and must be used together with other PPE protecting the uncovered parts of the body.

Utilization

In case of exposition to biological agents it must be used just once and it cannot be washed or be disinfected to use it again. If the PPE has been contaminated by any kind of liquid

and/or biological and chemical material, or a non visible contamination is supposed, change it.

ULTRASONIC AND
TAPED SEAMS

ELASTIC BAND

Tuta GVS2

GVS2 Coverall

TECHNICAL SPECIFICATIONS

(please read instruction for use for further information)

MEANINGFUL RESULTS CONCERNING THE BARRIER CHARACTERISTICS OF THE CLOTHING FOR THE PROTECTION TO BIOLOGICAL AGENTS (EN 14126/04)

Test pressure	0 Kpa	1.75 Kpa (0,25 psi)	3.50 Kpa (0,50 psi)	7.00 Kpa (1,00 psi)	14.00 Kpa (2.00 psi)
Sample 1	P	P	P	P	P
Sample 2	P	P	P	P	P
Sample 3	P	P	P	P	P

MEANINGFUL RESULTS CONCERNING THE BARRIER CHARACTERISTICS OF THE CLOTHING FOR THE PROTECTION AGAINST CHEMICAL AGENTS

Chemical product	Sample	Thickness (mm)	Passing through time (min)
Methotrexate	Average	0,27	> 480
5-Fluorouracil	Average	0,30	> 480
Vincristine Sulphate	Average	0,29	> 480
Daunorubicin Chlorhydrate	Average	0,29	> 480
Doxorubicin Chlorhydrate	Average	0,29	> 480
Hidroxite Sodium	Average	0,29	> 480
Sulphuric acid (96%)	Average	0,30	> 480
Cyclophosphamide	Average	0,27	> 480