Equorinometro

Equorinometro cod.9119

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Aequorinometer

Aequorinometer cod.9119

The Aequorin and GFP (Green Fluorescent Protein) are two protein expressed by jellyfish *Aequorea Visctoria*.

Both this protein there are an excellent tool to study the Ca^{++} mediated intracellular signals and pathway.

In fact is possible to express this protein as a Ca^{++} probe in a variety of cell types. Also, using the recombinant DNA technology, is possible to construct aequorin probes directed to specific cellular compartments and intracellular organelles

The Instrument is composed by a superior body containing a low noise photomultiplier connected to an **amplifier-discriminator** side of a dark box kept at 4°C. The photomultiplier is controlled by an electronic component that allow to regulate its sensibility through the voltage. The photomultiplier is placed in close proximity (2-3 mm) of a **perfusing chamber** (13 mm diameter, 2 mm height, volume 265 μ l). The perfusing chamber is sealed on the top with a cover slip held in place with a thin layer of silicon. Cells are continuously perfused via a **peristaltic pump** with medium thermostated via a water bath (**GT-80/P**) at 37°C. During manipulation of the cell chamber, the photomultiplier is protected from light by a shutter. During aequorin experiments, the shutter is opened and the chamber with cells is placed in close proximity of the photomultiplier. **The output of the amplifier-discriminators captured by a photon-counting board, transferred to a computer and stored for further analysis. The signal is immediately displayed on the computer through the connection with a graphic recorder 05.80L Linseeds.**

Complete aequorinometer:

- Photomultiplier H7360-01
- Termostic group GT-80/P
- Photon counter C8855-01
- Peristaltic Pump
- 220/110 volts

Thermostatic Group GT-80/P